



Discontinued Models

We are fortunate to live in a time of constant technological advances, enabling us to improve devices to help us work safer and more efficiently.

The Sedation Stethoscope® - Piezo version 5.1 is the most current version. As of May 1, 2024, we are no longer able to service the following amplifier versions:



Blue Classic Wireless
v. 2.0
2009



Blue Piezo Wireless
v. 2.0
2009



Green Classic Wireless
v. 3.0
2014




Green Piezo Wireless
v. 3.0
2014



White label but NO Power Switch on side
v. 4.0
2018

Amplifier and Chestpiece Cable Compatibility

The Chestpiece Cable was redesigned in early 2022 to enhance durability and connectivity. The chestpiece jack also included a new raised silver connector. The newly designed Chestpiece Cable is only compatible with v. 5.1 amplifiers.



Current Version



Our current Piezo version 5.1 includes an Amplifier with a charging port on the left side and a power switch on the right side. The bottom of the amplifier has a chestpiece port on the left, a recessed pairing button in the center, and a headset port on the right. A Bluetooth® Headset is included with up 18 hrs. of use. The Chestpiece Cable has increased durability and connectivity.

The Sedation Stethoscope® - Piezo version 5.1 can be used for 18-20 hrs. when fully charged.



TRADE IN

TRADE UP

PROGRAM

As a loyal Sedation Stethoscope® user, you qualify for an exclusive trade-in discount.

If you have any of the discontinued amplifiers shown above that are no longer serviceable or are not currently working, you may return the entire Sedation Stethoscope® unit (amplifier, chestpiece cable, and headset) to Sedation Resource for \$150 off a new unit.



We will provide a postage paid envelope for your return and take care of the recycling.

Use code TRADEIN to order your new Sedation Stethoscope® today for \$575 (reg. price \$725), item # 21-100-01S.